

Just Say NO to H₂0

Dihydrogen Monoxide is out there... Protect yourself & your Organic Synthesis Experiments with DriSolv® anhydrous solvents in specially packaged 100 ml septa sealed bottles with this offer from EMD Millipore.

The Ingenious DriSolv® Bottle

Purseal® elastomeric Septa reseals itself after each syringe withdrawal. Offers superior protection from moisture. Chemically Inert. It's even resistant to THF.

Specially threaded cap plug, allows additional protection from moisture. Can also accept exclusive EMD Millipore withdrawal tubes for nitrogen purge.

Buy 3 Cases, Get 1 Free plus a

Support the Dihydrogen Monoxide Cause by participating in our special offer. Get an amazing low moisture product and a Tee shirt that both say NO to H2O.

How EMD Millipore and DriSolv anhydrous solvents can help.

Buy three cases (6 x100 ml/cs) of DriSolv® anhydrous solvent in convenient septa sealed 100 ml bottles and receive one additional case (6 x100 ml) plus a special limited edition tee shirt absolutely free.


DriSolv[®] anhydrous solvents from EMD Millipore have water levels as low as 30 ppm. Our unique

Purseal[®] Septa reseals each time an aliquot of solvent is removed via syringe. In addition, only DriSolv[®] anhydrous solvent offers a secondary means of protecting your product with a unique threaded plug that screws into the over cap. This same over cap can also accept exclusive EMD Millipore accessories that allow solvent withdrawal under a nitrogen purge.

Redemption Instructions

Participating Products

Solvent	Size	Item	H ₂ 0 (ppm)
Acetonitrile	100ml	AX0143-7	<= 50
Alcohol Rgt	100ml	AX0447-7	<= 50
Benzene	100ml	BX0214-7	<= 30
Chloroform (alkene stb.)	100ml	CX1057-7	<= 50
Chloroform (ethanol stb.)	100ml	CX1056-7	<= 50
1,2 Dichloroethane	100ml	DX0794-7	<= 50
Dichoromethane	100ml	DX0834-7	<= 50
p-Dioxane	100ml	DX2092-7	<= 50
N,N Dimethylformamide	100ml	DX1727-7	<= 50
Ethyl Ether	100ml	EX0194-7	<= 50
Ethyl Acetate	100ml	EX0237-7	<= 50
n-Hexane	100ml	HX0304-7	<= 50
Isopropyl Alcohol	100ml	PX1827-7	<= 50
Methanol	100ml	MX0472-7	<= 50
Methyl-t-Butyl-Ether	100ml	MX0824-7	<= 50
1-Methyl 2-Pyrrolidone	100ml	MX1396-7	<= 50
Methylsulfoxide	100ml	MX1457-7	<= 50
Pyridine	100ml	PX2012-7	<= 50
Tetrahydrofuran (unstab.)	100ml	TX0277-7	<= 50
Tetrahydrofuran (BHT stab.)	100ml	TX0284-7	<= 50
Toluene	100ml	TX0732-7	<= 50



How to receive your free solvents & tee shirt.

Visit our website at www.emdchemicals. com/drisolvpromo.

Submit your claim on line along with a copy of your invoice and reference our promo code: DRISOLV4ME.

-Please allow 4-6 weeks for delivery.

Terms & Conditions

- Promotion valid from 2-1-12 thru 7-31-12.
- Only one redemption per customer
- Purchase on one order.
- Full cases only.

For more information, or to order product Contact us at:

Phone : 1-800-645-5476 Email: nacustomerservice@merckgroup.com

EMD Millipore 290 Concord Road Billerica, MA 01821 www.emdchemicals.com



www.emdchemicals.com