

**University of Michigan Biomedical Research Store Locations:
2297 Life Sciences Institute, 734.615.5701
2580 MSRBII, 734.615.5447**

Products stocked from PIERCE, a part of ThermoFisher Scientific

Cat #	Description	Unit	Storage
Chemiluminescent Detection Substrates			
PI 32209	Pierce ECL Western Blotting Substrate	250ml Kit	4 ⁰ C
PI 32106	Pierce ECL Western Blotting Substrate	500ml Kit	4 ⁰ C
PI 34080	SuperSignal ^R West Pico Substrate	500 ml	4 ⁰ C
PIA 34075	SuperSignal ^R West Dura, Extended Duration	100 ml	4 ⁰ C
PI 34095	SuperSignal ^R West Femto, Maximum Sensitivity	100 ml	4 ⁰ C
PI 34094	SuperSignal ^R West Femto, Trial Kit	30 ml	4 ⁰ C
Lysis Buffer			
PI 89900	RIPA Buffer	100 ml	4 ⁰ C
Protease & Phosphatase Inhibitors			
PI 78410	Halt TM Protease Inhibitor (2ml, 100x + 2ml EDTA)	Kit	4 ⁰ C
PI 78420	Halt TM Phosphatase Inhibitor (1ml, 100x)	1 ml	4 ⁰ C
Protein Prestained M.W. Markers			
PI 26691	TriChromRanger TM Prestained M.W. Markers, 48-well plate	1 plate	r.t.
PI 26681	BlueRanger ^R Prestained M.W. Markers, 48-well plate	1 plate	r.t.
In-Gel Stains			
PI 24590	GelCode Blue Stain Reagent	500 ml	4 ⁰ C
PI 24615	Imperial Protein Stain	1-Liter	r.t.
Common Blocking Buffers			
PI 37535	SuperBlock (TBS) Blocking Buffer	1-Liter	4 ⁰ C
PI 37536	SuperBlock T20 (TBS) Blocking, w/ Tween-20	1-Liter	4 ⁰ C
EN 37543	StartingBlock T20 (TBS) Blocking, with Tween-20	1-Liter	4 ⁰ C
Convenient Wash Buffer Packs			
PI 28372	BupH TM Phosphate Buffered Saline (each pack = 500ml)	40 packs	r.t.
PI 28379	BupH TM Tris Buffered Saline (each pack = 500ml)	10 packs	r.t.
Stripping Buffer			
PI 21059	Restore West Blot Stripping Buffer	500 ml	r.t. or 4 ⁰ C
Protein Assay Reagents			
PI 23200	Comassie Assay Kit, w/ Albumin Stds	950 ml	r.t. and 4 ⁰ C
PI 23236	Comassie Plus Assay Kit, w/ Albumin Stds	950 ml	r.t. and 4 ⁰ C
PI 23238	Comassie Plus Assay Kit	300 ml	4 ⁰ C
PI 23227	BCA Protein Assay Kit (250 std tube assays)	Kit	r.t.

Special Orders for all PIERCE Products are encouraged from the Biomedical Research Stores

If you have any questions, please contact Gary DeBusscher, your Key Account Manager,
(734) 657-8669, or gary.debusscher@thermofisher.com

